

**Is with great pleasure that we'd like to announce the release of CharruaPACS advanced 6.29**

In this new version we have made the following upgrades from our last version:

1. Added the Cornerstone DICOM HTML5 viewer for WEB browsers.
2. Have added encryption option for both DICOM TLS and WEB HTTPS.
3. Added preliminary support for DICOMWeb Rest API: WADO-RS, STOW-RS and QIDO-RS.
4. Have an Active Directory authentication option.
5. Have corrected several bugs founded on the previous version.

We would like to have your feedback on this new version, in case of updating an old installation, please contact us for the correct upgrade procedure.

[www.charruasoft.com](http://www.charruasoft.com)